

# SHAKER ORBITAL

axyos

## SO4136

<b>Shaker Type</b>	Reciprocating
<b>Tray dimensions (W x D mm)</b>	414x360
<b>Bench Dimensions (W x D x H mm)</b>	454 x 360 x 181
<b>Orbit / Stroke (mm)</b>	Adjustable 20, 30, 40
<b>Speed rev per min / strokes per min</b>	Adjustable 20 - 200
<b>Micro Digital Control</b>	Included
<b>Digital Speed/Time Indication</b>	Included
<b>Countdown Timer</b>	Included
<b>Slow Speed Start Up</b>	Included
<b>Maintains Last Setting on Power Up</b>	Included
<b>Non slip tray</b>	Included

### Approx Number of Erlenmeyer Flasks

<b>2 litre</b>	4
<b>1 litre</b>	4
<b>500ml</b>	12
<b>250ml</b>	12
<b>150ml</b>	24
<b>125ml</b>	24
<b>50ml</b>	42
<b>Max Speed with non slip tray(flasks)</b>	140 rpm

<b>Electrical</b>	220-240v 50-60Hz single phase (alternative power ranges available)
-------------------	--

### Options Description

tu	Adjustable rack to suit a range of flasks and bottles so & sr series small
tt	Tray to accommodate various tulip clips so & sr series small
tu	Adjustable rack to suit a range of flasks and bottles so& sr series large
tt	Tray to accommodate various tulip clips so & sr series large
ts	Increase std top tray size of up to 100mm on each dimension to suit special requirements (pls advise new size)
t2.?	(2) tier (Dual level) non slip rack (advise height for 1st level)
x50.?	Tulip clip to hold 50ml flask eg: x50.5= 50ml x 5 clips
x100.?	Tulip clip to hold 100ml flask eg: x100.5 = 100ml x 5 clips
x250.?	Tulip clip to hold 250ml flask eg: x250.5 = 250ml x 5 clips
x500.?	Tulip clip to hold 500ml flask eg: x500.5 = 500ml x 5 clips
x1.?	Tulip clip to hold 1litre flask eg: x1.5= 1000ml x 5 clips
x2.?	Tulip clip to hold 2litre flask eg: x2.5 = 2000ml x 5 clips
x?.?	Special size flask clips up to 6 litre



Quality Endorsed Company  
300 8881 111 1221  
Brisbane, Australia

